

Device, operating system, and browser compatibility

The Change Healthcare Stratus Imaging Share system can be used on a device meeting the following minimum requirements. Mobile devices support use of the Stratus Imaging Share viewer only, and they are not supported for diagnostic image review.

$F_{\mathfrak{P}},\, C_{\psi^{+},\, \Upsilon^{+}},\, H_{\psi^{+},\, \gamma^{+}\psi^{+},\, \gamma^{-}\psi^{+}},\, I_{\psi^{+},\, \Upsilon^{+},\, \Upsilon^{+}\psi^{+},\, \gamma^{+}\psi^{+}},\, I_{\psi^{+},\, \Upsilon^{+}\psi^{+},\, \gamma^{+}\psi^{+}},\, I_{\psi^{+},\, \Upsilon^{+}\psi^{+}},\, I_{\psi^{+},\, \Upsilon^{+}$

D _v /c v	O	M	B	N
Microsoft® Windows® PC	Microsoft® Windows® 7 or higher	8 GB RAM or greater	Current version of Google™ Chrome	Greater than 5 Mbps
Apple Mac™	OS X 10.14 (Windows Mojave) or higher	8 GB RAM or greater	Current version of Google™ Chrome	Greater than 5 Mbps

D _v /c v	Q	М	В	N
Apple™ mobile device	iOS 10.0 or higher	1 GB RAM or greater	Current version of Google™ Chrome or the current version of Apple Safari	Minimum: 5 Mbps or greater Recommended: 50 Mbps or greater
Android™ mobile device	Android 9.0 (Pie) or higher	1 GB RAM or greater	Current version of Google™ Chrome	Minimum: 5 Mbps or greater Recommended: 50 Mbps or greater

Optum

